

BC-5D

HEMATOLOGY CONTROLS

CONTROL

ASSAY VALUES AND EXPECTED RANGES

LOT

BC2209B**2022-11-10**

Instrument	Parameter	Low		Normal		High		+
		LOT	BC2209BL	LOT	BC2209BN	LOT	BC2209BH	
BC-5800, BC-5600	WBC ×10 ⁹ /L	3.47	± 0.50	8.06	± 1.00	18.11	± 2.50	
QC Mode	Neu# ×10 ⁹ /L	1.64	± 0.32	4.49	± 0.73	11.52	± 1.63	
	Lym# ×10 ⁹ /L	1.33	± 0.32	2.33	± 0.65	3.62	± 1.45	
	Mon# ×10 ⁹ /L	0.29	± 0.21	0.57	± 0.49	1.16	± 0.91	
	Eos# ×10 ⁹ /L	0.17	± 0.14	0.59	± 0.49	1.63	± 1.27	
	Bas# ×10 ⁹ /L	0.03	± 0.03	0.08	± 0.08	0.18	± 0.18	
	Neu%	47.3	± 9.0	55.7	± 9.0	63.6	± 9.0	
	Lym%	38.3	± 9.0	28.9	± 8.0	20.0	± 8.0	
	Mon%	8.5	± 6.0	7.1	± 6.0	6.4	± 5.0	
	Eos%	4.9	± 4.0	7.3	± 6.0	9.0	± 7.0	
	Bas%	1.0	± 1.0	1.0	± 1.0	1.0	± 1.0	
	RBC ×10 ¹² /L	2.14	± 0.18	4.24	± 0.24	5.06	± 0.30	
	HGB g/L	58	± 4	133	± 6	171	± 8	
	HCT %	17.7	± 1.5	40.1	± 2.0	52.4	± 2.4	
	MCV fL	82.8	± 5.0	94.5	± 5.0	103.5	± 5.0	
	MCH pg	27.1	± 2.5	31.4	± 2.5	33.8	± 2.5	
	MCHC g/L	327	± 30	332	± 30	327	± 30	
	RDW-CV %	16.0	± 3.0	13.9	± 3.0	13.7	± 3.0	
	RDW-SD fL	47.2	± 10.0	48.6	± 10.0	53.1	± 10.0	
	PLT ×10 ⁹ /L	48	± 20	252	± 40	505	± 60	
	MPV fL	10.1	± 3.0	9.9	± 3.0	9.3	± 3.0	
	PCT %*	0.048	± 0.048	0.249	± 0.100	0.470	± 0.200	
	PDW*	17.0	± 3.0	16.2	± 3.0	15.7	± 3.0	
	P-LCC ×10 ⁹ /L	19	± 15	91	± 25	140	± 35	
	P-LCR %	39.8	± 10.0	36.2	± 10.0	27.8	± 10.0	
BC-5390	WBC ×10 ⁹ /L	3.25	± 0.50	7.85	± 1.00	17.45	± 2.50	
QC Mode	Neu# ×10 ⁹ /L	1.63	± 0.30	4.51	± 0.71	11.34	± 1.57	
	Lym# ×10 ⁹ /L	1.20	± 0.30	2.20	± 0.71	3.32	± 1.23	
	Mon# ×10 ⁹ /L	0.24	± 0.20	0.43	± 0.32	0.96	± 0.70	
	Eos# ×10 ⁹ /L	0.18	± 0.14	0.71	± 0.56	1.83	± 1.40	
	Bas# ×10 ⁹ /L	0.82	± 0.33	2.14	± 0.79	5.36	± 1.75	
	Neu%	50.0	± 9.0	57.5	± 9.0	65.0	± 9.0	
	Lym%	37.0	± 9.0	28.0	± 9.0	19.0	± 7.0	
	Mon%	7.5	± 6.0	5.5	± 4.0	5.5	± 4.0	
	Eos%	5.5	± 4.0	9.0	± 7.0	10.5	± 8.0	
	Bas%	25.3	± 10.0	27.3	± 10.0	30.7	± 10.0	
	RBC ×10 ¹² /L	2.07	± 0.18	4.18	± 0.24	5.02	± 0.30	
	HGB g/L	53	± 4	122	± 6	159	± 8	
	HCT %	17.0	± 1.5	38.5	± 2.0	50.7	± 2.4	
	MCV fL	82.0	± 5.0	92.0	± 5.0	101.0	± 5.0	
	MCH pg	25.6	± 2.5	29.2	± 2.5	31.7	± 2.5	
	MCHC g/L	312	± 30	317	± 30	314	± 30	
	RDW-CV %	15.0	± 3.0	13.5	± 3.0	13.5	± 3.0	
	RDW-SD fL	47.5	± 8.0	48.0	± 8.0	51.5	± 8.0	
	PLT ×10 ⁹ /L	46	± 20	248	± 40	502	± 60	
	MPV fL	12.9	± 3.0	12.6	± 3.0	11.9	± 3.0	

* For Research Use Only

Before using, refer to the instruction sheet for mixing directions.

All brands and products are trademarks or registered trademarks of their respective companies.

BC-5D

HEMATOLOGY CONTROLS

CONTROL

ASSAY VALUES AND EXPECTED RANGES

LOT


BC2209B
2022-11-10

Instrument	Parameter	Low		Normal		High		++
		LOT	BC2209BL	LOT	BC2209BN	LOT	BC2209BH	
BC-5390 CRP	WBC ×10 ⁹ /L		3.26 ± 0.50		7.82 ± 1.00		17.39 ± 2.50	
QC Mode	Neu# ×10 ⁹ /L		1.60 ± 0.30		4.59 ± 0.71		11.20 ± 1.57	
	Lym# ×10 ⁹ /L		1.25 ± 0.30		2.15 ± 0.63		3.32 ± 1.40	
	Mon# ×10 ⁹ /L		0.22 ± 0.17		0.41 ± 0.40		1.04 ± 0.87	
	Eos# ×10 ⁹ /L		0.19 ± 0.17		0.67 ± 0.47		1.83 ± 1.40	
	Bas# ×10 ⁹ /L		0.83 ± 0.33		2.13 ± 0.78		5.34 ± 1.74	
	Neu%		49.2 ± 9.0		58.7 ± 9.0		64.4 ± 9.0	
	Lym%		38.4 ± 9.0		27.5 ± 8.0		19.1 ± 8.0	
	Mon%		6.6 ± 5.0		5.2 ± 5.0		6.0 ± 5.0	
	Eos%		5.8 ± 5.0		8.6 ± 6.0		10.5 ± 8.0	
	Bas%		25.6 ± 10.0		27.3 ± 10.0		30.7 ± 10.0	
	RBC ×10 ¹² /L		2.09 ± 0.18		4.18 ± 0.24		5.00 ± 0.30	
	HGB g/L		54 ± 4		124 ± 6		161 ± 8	
	HCT %		16.8 ± 1.5		38.9 ± 2.0		50.9 ± 2.4	
	MCV fL		80.4 ± 5.0		93.0 ± 5.0		101.8 ± 5.0	
	MCH pg		25.8 ± 2.5		29.7 ± 2.5		32.2 ± 2.5	
	MCHC g/L		321 ± 30		319 ± 30		316 ± 30	
	RDW-CV %		15.8 ± 3.0		14.2 ± 3.0		14.0 ± 3.0	
	RDW-SD fL		45.1 ± 8.0		46.2 ± 8.0		50.5 ± 8.0	
	PLT ×10 ⁹ /L		43 ± 20		243 ± 40		494 ± 60	
	MPV fL		10.6 ± 3.0		10.4 ± 3.0		9.6 ± 3.0	
	PCT %*		0.046 ± 0.046		0.253 ± 0.100		0.474 ± 0.200	
	PDW*		16.3 ± 3.0		16.1 ± 3.0		16.0 ± 3.0	
	P-LCC ×10 ⁹ /L		14 ± 14		73 ± 25		116 ± 35	
	P-LCR %		33.7 ± 10.0		30.2 ± 10.0		23.4 ± 10.0	
BC-5300, BC-5100	WBC ×10 ⁹ /L		3.20 ± 0.50		7.68 ± 1.00		17.40 ± 2.50	
BC-5380, BC-5180	Neu# ×10 ⁹ /L		1.63 ± 0.29		4.55 ± 0.70		11.54 ± 1.57	
QC Mode	Lym# ×10 ⁹ /L		1.21 ± 0.29		2.17 ± 0.70		3.31 ± 1.40	
(Software version lower than 1.24.00.16860)	Mon# ×10 ⁹ /L		0.17 ± 0.17		0.31 ± 0.24		0.70 ± 0.70	
	Eos# ×10 ⁹ /L		0.19 ± 0.16		0.65 ± 0.46		1.86 ± 1.40	
	Bas# ×10 ⁹ /L		1.83 ± 0.32		5.11 ± 0.77		13.57 ± 1.74	
	Neu%		51.0 ± 9.0		59.3 ± 9.0		66.3 ± 9.0	
	Lym%		37.8 ± 9.0		28.2 ± 9.0		19.0 ± 8.0	
	Mon%		5.2 ± 5.0		4.0 ± 3.0		4.0 ± 4.0	
	Eos%		6.0 ± 5.0		8.5 ± 6.0		10.7 ± 8.0	
	Bas%		57.3 ± 10.0		66.5 ± 10.0		78.0 ± 10.0	
	RBC ×10 ¹² /L		2.08 ± 0.18		4.17 ± 0.24		5.00 ± 0.30	
	HGB g/L		54 ± 4		124 ± 6		159 ± 8	
	HCT %		17.5 ± 1.5		40.8 ± 2.0		52.5 ± 2.4	
	MCV fL		84.3 ± 5.0		97.8 ± 5.0		105.0 ± 5.0	
	MCH pg		26.0 ± 2.5		29.7 ± 2.5		31.8 ± 2.5	
	MCHC g/L		308 ± 30		304 ± 30		303 ± 30	
	RDW-CV %		15.8 ± 3.0		13.5 ± 3.0		13.8 ± 3.0	
	RDW-SD fL		58.5 ± 8.0		59.0 ± 8.0		63.3 ± 8.0	
	PLT ×10 ⁹ /L		42 ± 20		239 ± 40		485 ± 60	
	MPV fL		10.4 ± 3.0		10.2 ± 3.0		9.4 ± 3.0	
	PCT %*		0.046 ± 0.046		0.244 ± 0.100		0.456 ± 0.200	
	PDW*		16.1 ± 3.0		16.0 ± 3.0		15.9 ± 3.0	

* For Research Use Only

Before using, refer to the instruction sheet for mixing directions.

All brands and products are trademarks or registered trademarks of their respective companies.

BC-5D

HEMATOLOGY CONTROLS

CONTROL

ASSAY VALUES AND EXPECTED RANGES

LOT

BC2209B**2022-11-10**

Instrument	Parameter	Low		Normal		High		+++
		LOT	BC2209BL	LOT	BC2209BN	LOT	BC2209BH	
BC-5300, BC-5100	WBC ×10 ⁹ /L		3.20 ± 0.50		7.57 ± 1.00		17.18 ± 2.50	
BC-5380, BC-5180	Neu# ×10 ⁹ /L		1.59 ± 0.29		4.53 ± 0.69		11.22 ± 1.55	
QC Mode	Lym# ×10 ⁹ /L		1.21 ± 0.29		2.10 ± 0.61		3.33 ± 1.38	
(Software version 1.24.00.16860 or higher)	Mon# ×10 ⁹ /L		0.18 ± 0.16		0.30 ± 0.24		0.76 ± 0.70	
	Eos# ×10 ⁹ /L		0.21 ± 0.19		0.64 ± 0.53		1.87 ± 1.38	
	Bas# ×10 ⁹ /L		1.87 ± 0.32		5.11 ± 0.76		13.54 ± 1.72	
	Neu%		49.7 ± 9.0		59.9 ± 9.0		65.3 ± 9.0	
	Lym%		37.9 ± 9.0		27.7 ± 8.0		19.4 ± 8.0	
	Mon%		5.7 ± 5.0		3.9 ± 3.0		4.4 ± 4.0	
	Eos%		6.7 ± 6.0		8.5 ± 7.0		10.9 ± 8.0	
	Bas%		58.5 ± 10.0		67.5 ± 10.0		78.8 ± 10.0	
	RBC ×10 ¹² /L		2.10 ± 0.18		4.17 ± 0.24		4.98 ± 0.30	
	HGB g/L		54 ± 4		124 ± 6		160 ± 8	
	HCT %		17.4 ± 1.5		39.4 ± 2.0		51.5 ± 2.4	
	MCV fL		82.9 ± 5.0		94.6 ± 5.0		103.5 ± 5.0	
	MCH pg		25.7 ± 2.5		29.7 ± 2.5		32.1 ± 2.5	
	MCHC g/L		310 ± 30		314 ± 30		310 ± 30	
	RDW-CV %		15.6 ± 3.0		14.0 ± 3.0		13.8 ± 3.0	
	RDW-SD fL		54.8 ± 8.0		55.5 ± 8.0		60.7 ± 8.0	
	PLT ×10 ⁹ /L		41 ± 20		245 ± 40		492 ± 60	
	MPV fL		10.0 ± 3.0		10.0 ± 3.0		9.2 ± 3.0	
	PCT %*		0.041 ± 0.041		0.245 ± 0.100		0.453 ± 0.200	
	PDW*		16.1 ± 3.0		16.1 ± 3.0		16.0 ± 3.0	
BC-5000, BC-5150, BC-5120	WBC ×10 ⁹ /L		3.30 ± 0.50		7.67 ± 1.00		17.24 ± 2.50	
BC-5130, BC-5140, BC-5000VET	Neu# ×10 ⁹ /L		1.54 ± 0.40		4.31 ± 0.92		10.86 ± 2.07	
QC Mode	Lym# ×10 ⁹ /L		1.24 ± 0.30		2.09 ± 0.62		3.10 ± 1.21	
	Mon# ×10 ⁹ /L		0.32 ± 0.32		0.58 ± 0.58		1.22 ± 1.22	
	Eos# ×10 ⁹ /L		0.17 ± 0.17		0.58 ± 0.58		1.74 ± 1.74	
	Bas# ×10 ⁹ /L		0.03 ± 0.03		0.11 ± 0.11		0.31 ± 0.31	
	Neu%		46.6 ± 12.0		56.2 ± 12.0		63.0 ± 12.0	
	Lym%		37.6 ± 9.0		27.3 ± 8.0		18.0 ± 7.0	
	Mon%		9.7 ± 9.7		7.5 ± 7.5		7.1 ± 7.1	
	Eos%		5.1 ± 5.1		7.6 ± 7.6		10.1 ± 10.1	
	Bas%		1.0 ± 1.0		1.4 ± 1.4		1.8 ± 1.8	
	RBC ×10 ¹² /L		2.09 ± 0.18		4.21 ± 0.24		5.06 ± 0.30	
	HGB g/L		54 ± 4		125 ± 6		164 ± 8	
	HCT %		17.4 ± 1.5		39.6 ± 2.0		51.2 ± 2.4	
	MCV fL		83.1 ± 5.0		94.1 ± 5.0		101.2 ± 5.0	
	MCH pg		25.8 ± 2.5		29.7 ± 2.5		32.4 ± 2.5	
	MCHC g/L		311 ± 30		316 ± 30		320 ± 30	
	RDW-CV %		18.8 ± 3.0		16.3 ± 3.0		16.2 ± 3.0	
	RDW-SD fL		56.8 ± 8.0		56.4 ± 8.0		60.6 ± 8.0	
	PLT ×10 ⁹ /L		44 ± 20		247 ± 40		505 ± 60	
	MPV fL		11.9 ± 3.0		12.2 ± 3.0		11.2 ± 3.0	
	PCT %*		0.052 ± 0.050		0.301 ± 0.100		0.566 ± 0.200	
	PDW*		15.6 ± 3.0		16.1 ± 3.0		16.2 ± 3.0	
	P-LCC ×10 ⁹ /L**		18 ± 15		101 ± 25		169 ± 35	
	P-LCR %**		41.1 ± 10.0		40.8 ± 10.0		33.4 ± 10.0	

* For Research Use Only

** These parameters are not provided on BC-5000/BC-5000 Vet analyzers

Before using, refer to the instruction sheet for mixing directions.

All brands and products are trademarks or registered trademarks of their respective companies.

BC-5D

HEMATOLOGY CONTROLS

CONTROL

ASSAY VALUES AND EXPECTED RANGES



BC2209B
2022-11-10

Instrument	Parameter	Low		Normal		High		++++
		LOT	BC2209BL	LOT	BC2209BN	LOT	BC2209BH	
BC-5300Vet, BC-5100Vet QC Mode	WBC ×10 ⁹ /L	3.20	± 0.50	7.68	± 1.00	17.40	± 2.50	
	Neu# ×10 ⁹ /L	1.63	± 0.29	4.55	± 0.70	11.54	± 1.57	
	Lym# ×10 ⁹ /L	1.21	± 0.29	2.17	± 0.70	3.31	± 1.40	
	Mon# ×10 ⁹ /L	0.17	± 0.17	0.31	± 0.24	0.70	± 0.70	
	Eos# ×10 ⁹ /L	0.19	± 0.16	0.65	± 0.46	1.86	± 1.40	
	Neu%	51.0	± 9.0	59.3	± 9.0	66.3	± 9.0	
	Lym%	37.8	± 9.0	28.2	± 9.0	19.0	± 8.0	
	Mon%	5.2	± 5.0	4.0	± 3.0	4.0	± 4.0	
	Eos%	6.0	± 5.0	8.5	± 6.0	10.7	± 8.0	
	RBC ×10 ¹² /L	2.08	± 0.18	4.17	± 0.24	5.00	± 0.30	
	HGB g/L	54	± 4	124	± 6	159	± 8	
	HCT %	17.5	± 1.5	40.8	± 2.0	52.5	± 2.4	
	MCV fL	84.3	± 5.0	97.8	± 5.0	105.0	± 5.0	
	MCH pg	26.0	± 2.5	29.7	± 2.5	31.8	± 2.5	
	MCHC g/L	308	± 30	304	± 30	303	± 30	
	RDW-CV %	15.8	± 3.0	13.5	± 3.0	13.8	± 3.0	
	RDW-SD fL	58.5	± 8.0	59.0	± 8.0	63.3	± 8.0	
	PLT ×10 ⁹ /L	42	± 20	239	± 40	485	± 60	
	MPV fL	10.4	± 3.0	10.2	± 3.0	9.4	± 3.0	
	PCT %*	0.046	± 0.046	0.244	± 0.100	0.456	± 0.200	
PDW*	16.1	± 3.0	16.0	± 3.0	15.9	± 3.0		

* For Research Use Only

Before using, refer to the instruction sheet for mixing directions.

All brands and products are trademarks or registered trademarks of their respective companies.



Mindray

Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

Mindray Building, Keji 12th Road South, Hi-tech Industrial Park, Nanshan, ShenZhen 518057, P.R.China

Tel: +86 755 81888998

Fax: +86 755 26582680

EC	REP
----	-----

Shanghai International Holding Corp. GmbH (Europe)

Eiffestraße 80 20537 Hamburg, Germany

Tel: 0049-40-2513175

Fax: 0049-40-255726